



Warsaw, 7th January 2020

Dear Valuable Customer,

We kindly inform that regarding the statistical analysis of results obtained in quality control tests of the HM22-001-1 CLEANER FOR MYTHIC 22 reagent, the specification for the pH parameter in the above product will be changed from 8.40 - 8.60 to 8.30 - 8.60. A lower 8.30 pH value does not affect the reagent functionality.

The change in the Certificate of Analysis will be implemented for reagent lots manufactured after January 1st, 2020.

No action required, that note is for informational purpose only.

In case of any questions please do not hesitate to contact your Area Sales Manager or my person directly at sonia.konicz@cormay.pl

Your sincerely

Sonia Konicz
Hematology Product Manager

NIO_2001_02